

INSTRUCTION FOR USE

Product Name: Zenit Brush – Pure Kolinsky Dental Composite Brush





Product Code:

- PRD.01.ZBR.0
- PRD.01.ZBR.000
- PRD.01.ZBR.1
- PRD.01.ZBR.SET.100

Product Description

Zenit Brush which is produced from pure Kolinsky enables you make perfect composite retouches.

Technical Specifications

	Product Code	Description	Ø	Length
1	PRD.01.ZBR.0	Glaze Brush #0 (Round)	1.2	7.5
				
2	PRD.01.ZBR.000	Glaze Brush #000 (Round)	1.2	5.5
				
3	PRD.01.ZBR.1	Opaque Brush #01 (Flat)	2.6	6.0
				
	PRD.01.ZBR.SET.100	Glaze Brush #0 (Round) + Glaze Brush #000 (Round) + Opaque Brush #01 (Flat)		
				

Instruction for Use

Indications for Dental Composite Brushes: Use the fine #1 brush for fine detailing and application of tints and opaques
Use the #2 and #3 brushes to apply bonding and silanating agents and for cement removal prior to polymerization.

Package Details

- PRD.01.ZBR.0 : 1x Glaze Brush #0 (Round)
- PRD.01.ZBR.000 : 1x Glaze Brush #000 (Round)
- PRD.01.ZBR.1 : 1x Opaque Brush #01 (Flat)
- PRD.01.ZBR.SET.100 : 1x Glaze Brush #0 (Round)
1x Glaze Brush #000 (Round)
1x Opaque Brush #01 (Flat)