# Low to High

## Arm holder

The table's vertical range of motion is 525-1,505 mm for the contour type and 715-1,275 mm for the pedestal type. The table is held by a flexible arm, which can be moved with a force 50% lighter than the previous models.



Handle suited for touch-panel operation
The handle is designed not only for an easy
grab-and-move motion, but also for providing a resting
position for the hand to reduce stress and ensure reliable
touch-panel operation.





#### Holder's position

The holder, designed to be independent from the table, can be freely adjusted with a three-axis adjustment mechanism for the position best suited for each dentist.



\*The installation image for left-handed specification.

# Release from weight

## Place holder

The dentist/dental assistant can pick up a handpiece and start treating patients right away without the need to re-hold the handpiece. Since the roller guide holds the weight of the hose, the dentist/dental assistant is freed from burden from the shoulder to the hand. They don't need to look back at the placeholder when returning the handpiece.

# nide oth rolling motion of the roller frees st/dental assistant from the weight se when he/she is drawing out a

#### Hose guide

The smooth rolling motion of the roller frees the dentist/dental assistant from the weight of the hose when he/she is drawing out a handpiece. This design helps him/her to administer treatment over a long period without feeling too much strain.





# Touch your style

## Interface

The interface has been designed to improve efficiency and reduce stress during treatment. When a handpiece is chosen and picked up, the screen displays only necessary information, helping the dentist concentrate on the treatment.



#### Motor

The setting screen for handpiece motors capable of storing four different settings allows intuitive operation. Rpm, torque and rotational direction can be chosen with one touch of the screen.



#### Scaler

From the handpiece setting screen, the treatment modes (periodontal, root canal and scaling) can be chosen with one touch of the screen.

# Improved cleaning performance

## Stay Clean

The hygienic condition inside the chair unit is maintained by water line flushing, so that patient-friendly treatment can be administered. As an environmentally friendly option, the unit comes with a separator or a separator with an amalgam collector, which helps reduce facilitating everyday maintenance.



Dismantling of the rod holder The end parts of the rod can easily be dismantled for cleaning.



Detachable holder

The dentist outlets and the assistant outlets can both be detached for cleaning.



Detachable light cover
The lens cover can be
removed with a single lever
action for easy cleaning.



Detachable bowl
The ceramic bowl can easily
be detached for cleaning
and can therefore always be
in a clean condition.



Water bottle Option



Cuspidor
The sliding cuspidor cover facilitates daily maintenance.



Flushing system Option
The integrated handpiece flushing system inhibits biofilm build up and bacterial growth in the handpieces hoses.



Vacuum Line Clean System Option

Daily use of the integrated suction cleaning system prevents build-up of debris in the vacuum line and helps to maintain constant suction.

## Comfort

## Stool

The stool goes well with Eurus. The seat and backrest can be adjusted for height and tilt.

(The assistant stool has a tilt function only.)

Meet the needs of various treatment styles
With the bucket-shaped seat and backrest, the stool is
fitted for treatment over a long period. With a low starting
height, the stool meets the needs of various treatment
styles.



## Color variations

## Leather Color

The colors come in four themes: warm, cool, earth and vitamin. Any choices from these leather colors should match the interior appointments of your clinic. The upholstery can be chosen from seamless, seam and luxury type.

#